

Name:

date:

rack #	6
1	D.I. water
2	Standard - 1.2
3	Sample - FS, smart water 1:10
4	Standard - 2.5
5	Sample - FS, Hewitt Tap 1:10
6	Sample - NY, 66 BF Station 1:10

rack #	5
1	D.I. Water
2	Standard .3
3	Sample - SH - HR 4:45 9/23 - 1:10
4	Standard 5
5	Sample - SH - HR 1:45 9/27 - 1:10
6	Sample - EB The Lake Central Park 1:10

rack #	3
1	D.I. water
2	Standard - 5
3	Sample NY 66 Tap Water 1:10
4	Standard - 1.2
5	Sample NY 66 Brita filter tap 1:10
6	Sample - RKS - AH tap 1:10

rack #	9
1	DI water
2	Standard - 5
3	Sample - EB Conservatory water central park 1:10
4	Standard - 1.2
5	Sample - EB - Harlem Meer central park 1:10
6	Sample - EB - The Pond central park 1:10

rack #	2 2
1	DI water
2	Standard - 2.5
3	RKS - tap water WH 1:10
4	TA - 1:40 NJ Tap H ₂ O
5	RKS - tap water Ms H 1:10
6	TA - fig H ₂ O (1:10)

rack #	10
1	DI Water
2	Standard - .3
3	Sample - KN - Corner Tap 1:10
4	Standard - 5
5	Sample - KN - 115 and Bway 1:10
6	Sample - AJC - EC faucet 1:10

rack #	8
1	DI water (blank)
2	Standard - 0.3
3	AL - croton reservoir 1:10
4	standard - 0.6
5	AL - Croton reservoir 1:50
6	AL - croton tap (ground) H ₂ O 1:10

rack #	12
1	DI Water
2	standard .6
3	Sample - AJC - puddle 1:10
4	CAL S12 Tap water 1:10
5	Sample - AJC - rain 1:10
6	Standard 2.5

rack #	7
1	DI Water (blank)
2	standard 5
3	SG puddle 1:50
4	Standard 2.5
5	SG pond 1:50
6	ACT NYCTAP 1:10

rack #	4
1	DI Water
2	Standard 2.5
3	Sample - KN - 115 and Bway 1:10
4	CAL Vapor distilled bottled H ₂ O 1:10
5	LI tap water - KD - 5:54 pm 9/26 1:10
6	Standard 0.3

↓
low
stds

rack #	1
1	DI WATER
2	Standard - 1.2
3	ACT CRANBURY TAP 1:10
4	Standard .6
5	ACT MAPLEWOOD TAP 1:10
6	ACT Princeton TAP 1:10

rack #	9
1	DI Water
2	Standard 1.2
3	Sample - KN - 115 and Bway 1:10
4	CAL 1055 artesian H ₂ O 1:10
5	KB Poland Spring 1:10
6	Sample - KN - 115 and Bway 1:10

Name:

date:

rack #	2
1	DI water
2	Standard 1.2ml
3	TA Fiji H ₂ O 1:10
4	AL Croton Reserve ~ 2dd off shore
5	NY Brita Filtered water
6	TA Fiji H ₂ O 1:10

rack #	
1	
2	
3	
4	
5	
6	

rack #	3
1	DI water
2	Standard 2.5ml
3	NY 616 Lobby filtered water bottle stat
4	RK's Jap water Alt 9/28/09
5	ES Smart H ₂ O
6	TA Fiji H ₂ O 1:10

rack #	
1	
2	
3	
4	
5	
6	

rack #	8
1	DI water
2	Standard 0.3
3	Std 0.6
4	Std 1.2
5	Std 2.5
6	Std 5

rack #	
1	
2	
3	
4	
5	
6	

rack #	
1	
2	
3	
4	
5	
6	

rack #	
1	
2	
3	
4	
5	
6	

rack #	
1	
2	
3	
4	
5	
6	

rack #	
1	
2	
3	
4	
5	
6	

rack #	
1	
2	
3	
4	
5	
6	

rack #	
1	
2	
3	
4	
5	
6	